MAYO CLINIC LABORATORIES

Test Definition: FBIUR

Bicarbonate, Urine

Overview

Method Name
Enzymatic Assay (EZA)

NY State Available
Yes

Specimen

Specimen Type
Urine

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Container/Tube: Plastic preservative-free urine container
Specimen Volume: 1 mL
Collection Instructions:
1. Collect 1 mL from a random urine collection in a preservative free plastic urine container.
2. Send refrigerated.

Specimen Minimum Volume
0.3 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	30 days	
	Ambient	48 hours	
	Frozen	30 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog.

Reference Values

Reporting limit determined each analysis.

None detected

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 11 days

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

82374

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBIUR	Bicarbonate, Urine	1964-6

Result ID	Test Result Name	Result LOINC® Value
Z3506	Bicarbonate	1964-6