MAYO CLINIC LABORATORIES

Test Definition: FONS

Western blot for anti-optic nerve autoantibodies in the serum

Overview

Special Instructions

- OHSU Requisition Form

Method Name

Western blot

NY State Available

Yes

Specimen

Specimen Type

Varies

Specimen Required

Complete and submit with specimen:
-Completed OHSU Ocular Requisition Form
-Clinical history
-Healthcare professional information (name and phone number)
NOTE: Without this information, testing cannot be completed.

Submit only 1 of the following specimens:

Specimen Type: Serum
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot serum into a plastic vial.
2. Label specimen as serum.
3. Send refrigerated.

Specimen Type: Plasma
Collection Container/Tube: Lavender top (EDTA)
Submission Container/Tube: Plastic vial

Specimen volume: 5 mL Plasma

Collection Instructions:

1. Centrifuge and aliquot plasma into a plastic vial.
2. Label specimen as plasma.
3. Send refrigerated.

Specimen Minimum Volume

3 mL

Reject Due To

Hemolysis	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated	7 days	

Clinical & Interpretive

Clinical Information

Refer to www.ohsu.edu/casey-eye-institute/ocular-immunology-lab

Reference Values

A final report will be provided.

Performance

PDF Report

Referral

Day(s) Performed

Batched

Report Available

16 to 35 days

Performing Laboratory Location

Ocular Immunology Laboratory OHSU

Test Definition: FONS

Western blot for anti-optic nerve autoantibodies in the serum

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

84182

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FONS	Anti-optic nerve autoantibodies, WB	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FONS	Anti-optic nerve autoantibodies, WB	Not Provided