MAYO CLINIC LABORATORIES

# Test Definition: FFTMV

Trichophyton Mentagrophytes (var interdigitale) IgE

## Overview

### Method Name
ImmunoCAP FEIA

### NY State Available
Yes

## Specimen

### Specimen Type
Serum

### Specimen Required
**Collection Container/Tube:**
**Preferred**: Red top
**Acceptable**: Serum gel
**Submission Container/Tube:** Plastic vial
**Specimen Volume:** 0.5 mL Serum
**Collection Instructions:**
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerate.

### Specimen Minimum Volume
See Specimen Required

### Reject Due To

| | |
|---|---|
| Lipemia | Reject |

### Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---|---|---|---|
| Serum | Refrigerated (preferred) | 28 days | |
| | Ambient | 28 days | |
| | Frozen | 28 days | |

## Clinical & Interpretive

# Test Definition: FFTMV

Trichophyton Mentagrophytes (var interdigitale) IgE

### Clinical Information

Refer to www.eurofins-viracor.com/test-menu

### Reference Values

<0.35 kU/L

### Interpretation

| Class | IgE (kU/L) | Comment |
|---|---|---|
| 0 | <0.10 | Negative |
| 0/1 | 0.10-0.34 | Equivocal/Borderline |
| 1 | 0.35-0.69 | Low Positive |
| 2 | 0.70-3.49 | Moderate Positive |
| 3 | 3.50-17.49 | High Positive |
| 4 | 17.50-49.99 | Very High Positive |
| 5 | 50.00-99.99 | Very High Positive |
| 6 | >99.99 | Very High Positive |

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

3 to 6 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

### Test Classification

# Test Definition: FFTMV

Trichophyton Mentagrophytes (var interdigitale) IgE

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

## CPT Code Information

86003

## LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---|---|---|
| FFTMV | Trichophyton ment var interdig IgE | Not Provided |

| Result ID | Test Result Name | Result LOINC® Value |
|---|---|---|
| Z5619 | Trichophyton ment var interdig IgE | 25801-2 |
| Z5620 | Class | 102906-5 |