MAYO CLINIC LABORATORIES

Test Definition: FCCA3

Cashew Component rAna o 3

Overview

Method Name
Phadia ImmunoCAP

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red or SST
Specimen Volume: 0.5 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume
0.5 mL

Reject Due To

Thawing:	Warm OK; Cold OK
Lipemia:	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values
<0.10 kU/L

Interpretation

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10 - 0.34	Equivocal/Borderline
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	High Positive
4	17.50 - 49.99	Very High Positive
5	50.00 - 99.99	Very High Positive
6	>99.99	Very High Positive

Clinical Reference

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
3 to 6 days

Performing Laboratory Location
Eurofins Viracor

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
86008

LOINC® Information

Test Definition: FCCA3

Cashew Component rAna o 3

Test ID	Test Order Name	Order LOINC® Value
FCCA3	Cashew Component rAna o 3	82539-8

Result ID	Test Result Name	Result LOINC® Value
Z5621	Cashew Component rAna o 3	82539-8
Z5622	Class	103115-2