MAYO CLINIC LABORATORIES

Test Definition: FASAB

Aspergillus Antibodies, Quantitative, DID

Overview

Method Name
Double Immunodiffusion (DID)

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube: Red top or serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Centrifuge immediately after coagulation (30 minutes) to prevent hemolysis.
2. Aliquot 1 mL of serum into a plastic vial.
3. Send frozen.

Specimen Minimum Volume
Serum: 0.5 mL

Reject Due To

Hemolysis	Reject
Lipemia	Reject
Gross bacterial contamination	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Ambient	14 days	
	Refrigerated	14 days	

Clinical & Interpretive

Clinical Information

Refer to www.labcorp.com/test-menu/

Reference Values

Aspergillus fumigatus
Neg: <1:1

Aspergillus flavus
Neg: <1:1

Aspergillus niger
Neg: <1:1

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

4 to 9 days

Performing Laboratory Location

LabCorp Burlington

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

86606 x 3

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FASAB	Aspergillus Ab, Qn, DID	Not Provided

Test Definition: FASAB

Aspergillus Antibodies, Quantitative, DID

Result ID	Test Result Name	Result LOINC® Value
Z5651	Aspergillus fumigatus	23821-2
Z5652	Aspergillus flavus	23820-4
Z5653	Aspergillus niger	10894-4