MAYO CLINIC LABORATORIES

Test Definition: FADPF

Adenosine Deaminase, Pleural Fluid

Overview

Useful For
Evaluation of tuberculous pleuritis

Method Name
Quantitative Spectrometry

NY State Available
Yes

Specimen

Specimen Type
Pleural Fluid

Specimen Required
Specimen Type: Pleural Fluid
Container/Tube: Plastic vial
Specimen Volume: 0.5 mL
Collection Instructions:
1.Collect pleural fluid in a leak proof container.
2. Centrifuge specimen at ambient temperature.
3. Transfer 0.5 mL into a plastic vial and freeze.
3. Ship frozen.
Note: Specimen must remain frozen until received at performing lab.

Specimen Minimum Volume
0.2 mL

Reject Due To

Whole blood	Reject
Bronchoalveolar lavage	Reject
Turbid specimen	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container

Pleural Fluid	Frozen (preferred)	30 days	
	Refrigerated	7 days	

Clinical & Interpretive

Reference Values
0-30 U/L

Performance

PDF Report
No

Day(s) Performed
Tuesday, Thursday, Sunday

Report Available
3 to 8 days

Performing Laboratory Location
ARUP Laboratories

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

CPT Code Information
84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FADPF	Adenosine Deaminase Pleural Fld	35704-6

MAYO CLINIC LABORATORIES

Test Definition: FADPF

Adenosine Deaminase, Pleural Fluid

Result ID	Test Result Name	Result LOINC® Value
Z5889	Adenosine Deaminase Pleural Fld	35704-6