MAYO CLINIC LABORATORIES

Test Definition: FIXCF

Immunofixation, CSF

Overview

Method Name
Immunofixation Electrophoresis

NY State Available
Yes

Specimen

Specimen Type
CSF

Specimen Required
Specimen Type: CSF
Container/Tube: sterile screw cap container
Specimen Volume: 6 mL
Collection Instructions: Collect 6 mL of spinal fluid (CSF) in a sterile screw cap container. Ship frozen.

Specimen Minimum Volume
4.5 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Ambient	6 days	
	Refrigerated	10 days	

Clinical & Interpretive

Reference Values
No abnormal bands are present on immunofixation.

Performance

Test Definition: FIXCF

Immunofixation, CSF

PDF Report
Referral

Day(s) Performed
Sunday through Friday

Report Available
5 to 9 days

Performing Laboratory Location
Quest Diagnostics Nichols Institute

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
86335

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIXCF	Immunofixation, CSF	13174-8

Result ID	Test Result Name	Result LOINC® Value
FIXCF	Immunofixation, CSF	13174-8