MAYO CLINIC LABORATORIES

Test Definition: FCDU6

Cocaine and Metabolites, Umbilical Cord Tissue (Non-Orderable)

Overview

Useful For

Only orderable as a profile. For more information see FCDSU / Comprehensive Drug Screen, Umbilical Cord Tissue.

Testing Algorithm

If appropriate, confirmation testing will be performed at no additional charge.

Method Name

Only orderable as a profile. For more information see FCDSU / Comprehensive Drug Screen, Umbilical Cord Tissue.

NY State Available

Yes

Specimen

Specimen Type

Tissue

Specimen Required

Only orderable as a profile. For more information see FCDSU / Comprehensive Drug Screen, Umbilical Cord Tissue.

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Tissue	Refrigerated (preferred)	7 days	
	Frozen	14 days	

Clinical & Interpretive

Clinical Information

Refer to www.nmslabs.com/test-catalog.

Reference Values

Only orderable as a profile. For more information see FCDSU / Comprehensive Drug Screen, Umbilical Cord Tissue.

Reporting limit determined each analysis.

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
3 to 7 days

Performing Laboratory Location
NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCDU6	Cocaine and Metabolites, UC	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6176	Benzoylecgonine	97222-4
Z6177	Cocaethylene	97232-3
Z6178	Cocaine	97234-9