MAYO CLINIC LABORATORIES

# Test Definition: FPROS

Prostaglandin D2 (PG D2), Random Urine

## Overview

### Method Name
Direct ELISA

### NY State Available
Yes

## Specimen

### Specimen Type
Urine

### Specimen Required
**Patient Preparation:** Patient **should not** be on aspirin, indomethacin, or anti-inflammatory medications, if possible, for at least 48 hours prior to collection of specimen.
**Submission Container/Tube**: Plastic vial
**Specimen Volume**: 10 mL
**Collection Instructions**:
1. Collect 10 mL random urine (24-hour urine collection is **not acceptable**)
2. No preservative.
2. Freeze immediately and send specimen frozen in a plastic vial.
3. Specimen must remain frozen prior to and during shipping.

### Specimen Minimum Volume
5 mL

### Reject Due To

| | |
|---|---|
| Gross hemolysis | Reject |
| Gross lipemia | Reject |
| Gross icterus | Reject |
| 24-hour urine collection | Reject |
| Ambient or refrigerated specimens | Reject |

### Specimen Stability Information

# Test Definition: FPROS

Prostaglandin D2 (PG D2), Random Urine

| Specimen Type | Temperature | Time | Special Container |
|---|---|---|---|
| Urine | Frozen | 60 days | |

## Clinical & Interpretive

### Clinical Information

Refer to www.interscienceinstitute.com/individual-assays/

### Reference Values

Up to 175 ng/g Creatinine

## Performance

### PDF Report

Referral

### Day(s) Performed

Monday through Friday

### Report Available

5 to 9 days

### Performing Laboratory Location

Inter Science Institute

## Fees & Codes

### Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

### Test Classification

This test has not been cleared or approved by the US Food and Drug Administration.
This test was developed and its performance characteristics determined by Inter Science Institute. Values obtained with different methods, laboratories, or kits cannot be used interchangeably with the results on this report. The results cannot be interpreted as absolute evidence of the presence or absence of malignant disease.

### CPT Code Information

84150

# Test Definition: FPROS

Prostaglandin D2 (PG D2), Random Urine

## LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---|---|---|
| FPROS | Prostaglandin D2 (PG D2), U | 12838-9 |

| Result ID | Test Result Name | Result LOINC® Value |
|---|---|---|
| FPROS | Prostaglandin D2 (PG D2), U | 12838-9 |