MAYO CLINIC LABORATORIES

Test Definition: FFTSI

Thyroid-stimulating Immunoglobulin (TSI)

Overview

Method Name
IMMULITE 2000 TSI

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red-top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL serum
Collection Instructions:
1. Within 2 hours of collection, centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume
Serum: 0.3 mL Note: This volume does not allow for repeat testing

Reject Due To

Gross lipemia	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen	14 days	

Clinical & Interpretive

Clinical Information

Refer to https://www.labcorp.com/test-menu/

Reference Values
0.00-0.55 IU/L

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
4 to 8 days

Performing Laboratory Location
LabCorp Burlington

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
84445

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFTSI	Thyroid Stim Immunoglobulin	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FFTSI	Thyroid Stim Immunoglobulin	30567-2