Test Definition: FNNAT

Neonatal Alloimmune Thrombocytopenia (NAIT)-Serial Monitoring Maternal only

Overview

Special Instructions

- Versiti Platelet Neutrophil Immunology Lab Form

Highlights

This assay includes: Testing for antibodies against HPA-1a/b, HPA-2a/b, HPA-3a/b, HPA-4a, HPA-5a/b, GPIIb/IIIa, GPIa/IIa, GPIb/IX. GPIV and Class I HLA.

Method Name

Platelet Antibody Bead Array (PABA) and Flow Cytometry

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Specimen Type: Serum
Collection Container/Tube: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 5 mL
Collection Instruction:
1.Centrifuge and aliquot 5 mL serum into a plastic vial.
2.Send refrigerated.

Forms

Platelet Neutrophil Immunology Laboratory **is required**. Complete and send with the specimen.

Specimen Minimum Volume

3 mL

Reject Due To

Gross hemolysis	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values
Negative

Performance

PDF Report
Referral

Day(s) Performed
Varies

Report Available
9 to 11 days

Performing Laboratory Location
Versiti Blood Center of Wisconsin

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by Versiti Wisconsin, Inc. It has not been cleared or approved by the US Food and Drug Administration. This test is used for clinical purposes. It should not be regarded as investigational or for research. This laboratory is certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical laboratory testing.

CPT Code Information
86022

Test Definition: FNNAT

Neonatal Alloimmune Thrombocytopenia (NAIT)-Serial Monitoring Maternal only

86023 x 2

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FNNAT	Neonatal Alloimmune Thromb (NAIT)	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FNNAT	Neonatal Alloimmune Thromb (NAIT)	Not Provided