MAYO CLINIC LABORATORIES

Test Definition: FDISU

Diuretic Profile, Urine

Overview

Method Name
High Performance Liquid Chromatography with Ultraviolet Detection; (HPLC-UV)

NY State Available
Yes

Specimen

Specimen Type
Urine

Specimen Required
Specimen Type: Urine
Container/Tube: Plastic preservative-free urine container
Specimen Volume: 10 mL
Collection Instructions:
1. Collect 10 mL random urine without preservative.
2. Send specimen refrigerated.

Specimen Minimum Volume
1.2 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	14 days	

Clinical & Interpretive

Performance

Test Definition: FDISU

Diuretic Profile, Urine

PDF Report
No

Day(s) Performed
Thursday

Report Available
7 to 9 days

Performing Laboratory Location
Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
80377

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDISU	Diuretic Screen, Urine	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6390	Diuretic Screen	Not Provided
Z6391	Bumetanide	3409-0
Z6392	Chlorothiazide	9508-3
Z6393	Chlorthalidone	3478-5
Z6394	Furosemide	3660-8
Z6395	Hydrochlorothiazide	3676-4
Z6396	Hydroflumethiazide	40469-9
Z6397	Metolazone	12347-1