MAYO CLINIC LABORATORIES

Test Definition: FLXPN

Loxapine, Plasma

Overview

Method Name
High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Plasma EDTA

Specimen Required
Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Collection Container/Tube: Lavender or pink top (EDTA)
Submission Container/Tube: Plastic screw-capped vial
Specimen Volume: 1 mL plasma
Collection Instructions:
1. Promptly centrifuge and aliquot 1 mL plasma into a plastic screw-capped vial.
2. Send refrigerated.
Specimen Stability Information: Refrigerated (preferred) 30 days/Ambient 30 days/Frozen 15 months

Specimen Minimum Volume
Plasma: 0.5 mL

Reject Due To

Polymer gel separation tube (SST or PST)	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Refrigerated (preferred)	30 days	
	Ambient	30 days	
	Frozen		

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/tests

Reference Values
None Detected
Dose and route: peak plasma concentration
Average 29.5 mg (15-50 mg) IM: 16 ng/mL 2 hr post-dose
Average 29.5 mg (15-50 mg) oral: 16 ng/mL 1 hr post-dose
Single 50 mg oral: 33 ng/mL 1.6 hr post-dose
Chronic 150 mg oral: 30 ng/mL
Half-life: approximately 6.9 hr

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
9 to 11 days

Performing Laboratory Location
NMS Labs

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information
80342

Test Definition: FLXPN

Loxapine, Plasma

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FLXPN	Loxapine	3728-3

Result ID	Test Result Name	Result LOINC® Value
FLXPN	Loxapine	3728-3