MAYO CLINIC LABORATORIES

Test Definition: FAVI

Hypersensitivity Pneumonitis Avian Panel

Overview

Method Name
Immunodiffusion

NY State Available
No

Specimen

Specimen Type
Serum

Specimen Required
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 3 mL of serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume
Serum: 1 mL

Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Gross icterus	Reject
Empty container	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	30 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information

Refer to
www.mcw.edu/departments/pediatrics/divisions/asthma-allergy-and-clinical-immunology/patient-care/allergy-immunology-diagnostic-lab-center

Reference Values

The following antigens will have results reported:
Pigeon Sera
Pigeon DE
Cockatiel
Parakeet
Parrot

Results are reported as:
Negative
Weakly Positive
Positive
Strong Positive

This result must be correlated with patient's clinical response and should not solely be considered in the diagnosis.

Performance

PDF Report

No

Day(s) Performed

Thursday and Friday

Report Available

18 to 20 days

Performing Laboratory Location

Medical College of WI Allergy-Immunology MACC Fund Research Ctr.

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

86331x5

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FAVI	Avian Panel	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z0289	Pigeon DE	Not Provided
Z0285	Parakeet	Not Provided
Z0286	Cockatiel	Not Provided
Z0287	Parrot	Not Provided
Z0288	Pigeon Sera	Not Provided