MAYO CLINIC LABORATORIES

Test Definition: FMETX

Methsuximide (Celontin) as Desmethylmethsuximide

Overview

Method Name
Gas Chromatography/Flame Ionization Detection (GC-FID)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
Submit only 1 of the following specimens:

Plasma
Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable.** Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

Serum
Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable.** Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume
0.25 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

	Frozen	180 days	

Clinical & Interpretive

Reference Values
10.0 - 40.0 ug/mL

Methsuximide measured as desmethylmethsuximide.

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
3 to 7 days

Performing Laboratory Location
Medtox Laboratories, Inc.

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
80339

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FMETX	Methsuximide (Celontin)	3801-8

Result ID	Test Result Name	Result LOINC® Value

Z1107	Methsuximide	3801-8